
NanoCycler

Easy, compact, and absolutely low-cost battery charger-discharger by NANO**BASE**



<u>Higher</u> performance instrument for battery researchers

Although many advanced battery research starts with coin cells, most researchers are forced to use standard battery cyclers designed for testing large cells such as 18650 cells. Responding to researchers' demand of more cost-effective battery cyclers, Nanobase, Inc. has launched NanoCycler, the higher performance battery cycler for research of coin cells. NanoCycler is easy and affordable for any battery researchers in chemistry, chemical engineering, electric engineering.

Features

Smooth

Place your 2025, 2032 coin cells directly into NanoCycler

Simple

Just connect NanoCycler to PC via USB port for analysis

Compact

Using in constant temperature chambers? No problem!

Expandable

Channel addition available for more than 80 channels per PC

Specifications

Product Size

 $153 \times 270 \times 107 (mm)$

Electricity

- 88 ~ 264 VAC or 125 ~ 373 VDC
- 80 W (max)

Voltage

- Range: 0 ~ 5 V
- Accuracy: ±0.1% FS
- · Measurement resolution: 16 bit
- Programming resolution: 14 bit

Current

- Range: 3 manually selectable ranges (200 µA, 2 mA, 20 mA)
- Accuracy: ±0.1% FS
- Measurement resolution: 16 bit
- Programming resolution: 14 bit
- Custom current range available

Channels

- 8 independent channels per unit
- · Sockets for 2025, 2032 coin cells
- More than 80 channels per PC (depending on PC specifications and available USB ports)

Data Recording Rate

- 1 kHz (max)
- 0.001 s ~ 9999 s

Software

- Sequence editor: Step & loop sequence programming
- Channel monitor & control
- Channel summary
- Plotting function : general plot, cycle plot, real-time plot
- Data export in .csv format


Software Capabilities

Channel control and monitoring

Channe	lOsp.														
Ch#	Label	Contro	Real Time	Status	Cycle #	Step	T	I (mA)	Range	V (V)	Ch Q (mAh)	Dch Q (mAh)	Sequence file	Data file	
1		Start	V	STBY	1	3	0.00.00	-0.0005	2 mA	-0.0053	20.28270	0.00000	Li_1C_170uah	NanoCycler/Data/ch1a	
2		Stop	VS	DISC	38	1	0:04:28	0.0629	2 mA	0.0046	0.00000	0.00082	2200uF_10ms	NanoCycler/Data/ch2	
3		Stop	Ve	CHG	3	2	0:00:14	2.0011	2 mA	2.2271	0.00033	0.00083	2200uF_10ms	NanoCycler/Data/ch3	
4		Stop	V×	DISC	3	1	0.00.12	0.0809	2 mA	0.0061	0.00000 0.00082		2200uF_10ms	NanoCycler/Data/ch4	
		Stop	VS:	CHG			0:00:10		2 mA	5.0003			2200uF_10ms	NanoCycler/Data/ch5	
6		Start	V	STBY	10000	5	0:00:00	0.0000	2 mA	0.2677	0.00086	0.00086	2200uF_10ms	NanoCycler/Data/ch6	
7		Start	V	STBY	10000	5	0:00:00	-0.0006	2 mA	-0.0050	0.00087	0.00087	2200uF_10ms	NanoCycler/Data/ch7	
8		Start	V	STBY	10000	5	0:00:00	0.0007	2 mA	0.2620	0.00088	0.00088	2200uF_10ms	NanoCycler/Data/ch8	
9		Start	Ve	ERR	3	5	0:00:15	-00.103	20 mA	-0.0063	0.00292	0.00296	2200uF_1ms	NanoCycleriData/ch91	

- Channel Start / Stop control with password protection
- Displays Cycle No., Step No., Elapsed Time, Current, Voltage, Q, Sequence File, and Data File

Plotting



Real-Time Plot

Sequence Editing

/s	equence Edit	or												
	Туре		1	lunit		V (V)	Cut-off type		Cut-off cond		Cut-off Value	Goto Step	Loop	Param1
•	Standby	٠	1	C	٠	2.8	StepTime	٠	GreaterT	٠	10			
П	Discharge	•	1	C	٠	2,8	Voltage	*	LessThan	*	2,9	0	0	0
	Charge	٠	1	С		4.2	Current		LessThan	٠	0,05	0	0	
	Charge	•	1	С	٠	4.2	StepTime	*	GreaterT		1	0	0	
	Loop	٠					Current		LessThan		0,02	-3	50	

- Step types: Charge, Discharge, Standby, Loop
- Cut-off types: Voltage, Current, Step Time, Cycle Time, Capacity

Channel Summarizing



- Charge, Discharge, Standby/Error Status of all channels shown as icons
- Also displays Elapsed Time, Current, and Voltage